



Water Treatment for Hemodialysis

Shiva Seyrafian

Nephrologist

IKRC- Isfahan University of Medical Sciences

1404/2/9 – 29/ 4/2025


Water treatment system standards

Standards should be based upon:

- International Organization for Standardization (ISO) issued by
- American National Standards Institute (ANSI) with
- the Association for the Advancement of Medical Instrumentation (AAMI).
 - **Water quality requirements and the installation,**
 - **Disinfection,**
 - **Maintenance procedures for water treatment systems.**

2025 UpToDate, Inc. Topic 1896 Version 31.0

The waterborne contaminants and the safe levels in hemodialysis



water contaminants

Contaminant	ANSI/AAMI/ISO standard (mg/L)	Contaminant	ANSI/AAMI/ISO standard (mg/L)
Arsenic	0.005	Nitrate (N)	2
Antimony	0.006	Potassium	8
Barium	0.1	Selenium	0.09
Beryllium	0.004	Silver	0.005
Cadmium	0.001	Sodium	70
Chromium	0.014	Sulfate	100
Copper	0.1	Zinc	0.1

ANSI: American National Standards Institute; AAMI: American Association for the Advancement of Medical Instrumentation; ISO: International Organization for Standardization

water contaminants...

Contaminant	ANSI/AAMI/ISO standard (mg/L)
Fluoride	0.2
Lead	0.005
Mercury	0.0002
Aluminum	0.01
Calcium	2
Magnesium	4
Total chlorine (free chlorine and chloramine)	0.1

ANSI: American National Standards Institute; AAMI: American Association for the Advancement of Medical Instrumentation; ISO: International Organization for Standardization



آزمایشگاه شرکت آب و فاضلاب استان اصفهان
Isfahan Water and Wastewater Laboratory

اصفهان - خیابان هزارحریب - خیابان جابری حیان - شرکت آب و فاضلاب استان اصفهان
Isfahan water & wastewater Company, Jaber Ehn Hayyan St., Hezar
Jerib St., Isfahan



LQS-E-708-III

Rev: 01

تاریخ انتشار: ۱۳۹۴/۱۱/۰۱

گزارش آزمون

Test Report

E03047

1403/07/23

نتایج آزمون (Test Results)

حد مجاز Specification	استاندارد مرجع Reference Standard	عدم قطعیت اندازه گیری (%) Uncertainty	روش آزمون Test Method	واحد Unit	نتیجه آزمون Test results	نام آزمون (پارامتر) Test Name (Parameter)
6	ISIRI 1053	22.37	SM 3112 B	µg/L	Result<1.05	Hg
10	ISIRI 1053	18.09	SM 3113 B	µg/L	Result<2.06	Pb
50	ISIRI 1053	10.43	SM 3113 B	µg/L	Result<1.9	Cr
3	ISIRI 1053	21.46	SM 3113 B	µg/L	Result<0.23	Cd
10	ISIRI 1053	20.93	SM 3114 C	µg/L	Result<2.14	As

Clinical manifestations of exposure of hemodialysis patients to water contaminants



water treatment for hemodialysis

Clinical manifestations of water contaminants

Aluminum	Neurologic deterioration, cognitive changes, encephalopathy, bone disease, anemia, death
Chlorine	Hemolysis, hyperkalemia, shortness of breath, hypoxia, cardiac arrest
Chlorine dioxide	Likely to cause methemoglobinemia and hemolytic anemia
Chloramine, (monochloramine)	Hemolysis, anemia, methemoglobinemia, cyanosis, shortness of breath, EPO resistance
Copper	Nausea, anemia, hemolysis, chills, death
Zinc	Anemia, cardiovascular events

Clinical manifestations of water contaminants..

Fluoride	Bone disease, pruritus, chest pain, nausea, cardiac arrest (acute exposure to high concentrations)
Lead	Abdominal pain, muscle weakness, pruritus, peripheral neuropathy, cognitive changes
Nitrate	Nausea and vomiting, hypotension, methemoglobinemia
Bacteria, endotoxin	Hypotension, chills, fever, death
Sulfate	Nausea and vomiting, metabolic acidosis
Hydrogen peroxide	Cyanosis, methemoglobinemia
Silver-stabilized hydrogen peroxide	Cyanosis, methemoglobinemia

Clinical manifestations of water contaminants..

Sodium hypochlorite	Hemolysis, hyperkalemia, shortness of breath, hypoxia, cardiac arrest
Peracetic acid	Hemolysis
Cyanotoxins	Acute illness, death
Microcystins	Acute illness, Visual disturbances, Liver failure, death,
Calcium	Nausea, vomiting, weakness, hypertension (hard water syndrome)
Magnesium	Nausea and vomiting, muscle weakness
Potassium	Cardiovascular events
Sodium	Hypertension

Components and Monitoring of Water treatment systems



Components and Monitoring of Water Treatment Systems


- 1. Booster pump**
- 2. Ultrafilters, multimedia filters or sediment or depth filters**
- 3. Softeners**
- 4. Activated carbon beds or filters**
- 5. Reverse osmosis**
- 6. Deionization**
- 7. Distribution system components**
- 8. Disinfection**

Components and Monitoring of Water Treatment Systems

1. Booster pump


Booster pump: help regulate flow when the **water pressure** falls below predetermined levels.





Water Treatment Systems

1. Booster pump
- 2. Ultrafilters, particle filters, or sediment or multimedia filters**
3. Softeners
4. Activated carbon beds or filters
5. Reverse osmosis
6. Deionization
7. Distribution system components
8. Disinfection



2- Ultrafilters, particle filters or multimedia filters


- Filters should be **backwashed every day** when there are no patients connected to the hemodialysis machines.
- **Pre- and post-filter pressure gauges** should be fitted on all filters and pressure difference between the two gauges should be monitored.
- All types of filters become progressively **obstructed with debris**, resulting in an increase in **pressure difference** and/or a decrease in **filtrate flow** rate.

Ultrafilter for
endotoxin removal



Pressure gauge






20 PSI → 15 PSI

water treatment for hemodialysis

MULTI-MEDIA FILTERS (MMF) (Sediment filter)




Water Treatment Systems

1. Booster pump
2. Ultrafilters, particle filters, multimedia filters, and depth filters
- 3. Softeners**
4. Activated carbon beds or filters
5. Reverse osmosis
6. Deionization
7. Distribution system components
8. Disinfection

3- Softeners

- In which **calcium and magnesium are exchanged for sodium ions**.
When the water is "hard" due to high levels of calcium and magnesium.
- Softeners must be **regularly regenerated**.
- Regeneration on an automatic cycle **by exposing the resin to a strong brine solution**.
- To ensure that calcium and magnesium ions remove effectively, **at the end of the day**, hardness should be checked.


Water softener



Monitoring of Softeners

- Colorimetric test strips can also be used to measure water hardness on site.
- Softened water should have a hardness of 1 GPG (17 mg/L CaCO_3) or lower.
- The brine tank should be checked daily to ensure adequate quantities of sodium chloride.





Water Treatment Systems

1. Booster pump
2. Ultrafilters, particle filters, multimedia filters, and depth filters
3. Softeners
- 4. Activated carbon beds or filters**
5. Reverse osmosis
6. Deionization
7. Distribution system components
8. Disinfection

4- Activated carbon beds or filters

- Activated carbon bed is the principal method of removing **chlorine** and **monochloramine**, substances that:
 1. Can **degrade reverse osmosis** membranes, and
 2. Cause **hemolysis** and **resistance to erythropoiesis-stimulating agents** (ESAs).
- Carbon can remove a wide range of organic contaminants by adsorption.



GAC Carbon Media




Granular Activated Carbon
is Contained in Cannister
in Loose form.

CB Carbon Media



Powder Activated Carbon
is Compressed & Bonded
into a Solid Block form.






4- Activated carbon beds or filters..

Monitoring: to detect exhaustion of their capacity for removal of monochloramine and total chlorine.

- Total chlorine measurements can be made by using a "dip and read" test strip.
- Testing should be performed **at the beginning of the day** after allowing the system **to run for 15 minutes** and then at least **every three to four hours** while patients are on dialysis.

Water Treatment Systems

1. Booster pump
2. Ultrafilters, particle filters, multimedia filters, and depth filters
3. Softeners
4. Activated carbon beds or filters
- 5. Reverse osmosis**
6. Deionization
7. Distribution system components
8. Disinfection



5- Reverse osmosis

- Removes dissolved substances, such as metal ions, salts, and inorganic and organic chemicals, and microbial elements, including bacteria, endotoxins, and viruses.
- Organic contaminants with molecular weight >200 Da are effectively filtered out.
- Reverse osmosis reduces contaminant levels by using high pressure (200 to 250 pounds/in²) to force a portion of the feed water across a semipermeable membrane.
- The remaining feed water that does not go across the membrane carries the concentrated contaminants to the drain.

5- Reverse osmosis..



- Although there is potential for up to 75 percent of the feed water to be recovered as product water, **recovery rates of 35 to 50%** are more common
- In case of **deficient feed water**, a part of **the reject water backs to the feed water** to allow recovery of a greater quantity of product water.
- **Two-stage reverse osmosis** can achieve an efficiency of **90 percent** or more.

5- Reverse osmosis..

- The most common reverse osmosis membrane is made of polyamide. Polyamide membranes are sensitive to chlorine and monochloramine, and to peracetic acid (disinfectant) at a concentration of >1%.

Monitoring:

- Ensure continued removal of contaminants at the desired level
- Ensure production of sufficient product water
- Determine when the membranes should be cleaned or replaced



A



B

Water Treatment Systems

1. Booster pump
2. Ultrafilters, particle filters, multimedia filters, and depth filters
3. Softeners
4. Activated carbon beds or filters
5. Reverse osmosis
- 6. Deionization**
7. Distribution system components
8. Disinfection

6-Deionization


- Ion-exchange resins are used to **remove ionic contaminants** from water.
- Hydrogen ions are exchanged for cations (**cationic resins**) and hydroxyl ions are exchanged for anions (**anionic resins**).
- If **feed water high contaminants** and **RO** is insufficient, **deionization** can be **added** as secondary purification step.
- **Does not remove** uncharged contaminants, **bacteria, or endotoxins**.

6-Deionization..

- The resins provide a good environment for bacterial proliferation and often worsen the microbiologic quality of water passing through them.
- A bacteria- and endotoxin-retentive filter to prevent contamination.

Monitoring:

- Measure of the concentration of ions in the water continuously.



Water Treatment Systems

1. Booster pump
2. Ultrafilters, particle filters, multimedia filters, and depth filters
3. Softeners
4. Activated carbon beds or filters
5. Reverse osmosis
6. Deionization
- 7. Distribution system components**
8. Disinfection

7- Distribution system components

- Metals and chemicals leach from materials to water.
- **Surfaces** in contact with treated water must be **fabricated with inert materials**.
- **No brass, aluminum, or galvanized** metal parts can be used.
- Commonly used **materials** include cross-linked **polyethylene** (PEX), polyvinylchloride (**PVC**) and **Stainless Steel**.
- If **hot water** is used as the sterilant for the distribution system, **PVC should not be used** due to its degradation following repeated hot water exposure.



**Polyethylene (PEX)
Tube**



PVC pipes



Stainless steel pipe

7- Distribution system components..

- **Drinking water** supply source independent from other water supply sources on site.
- **The storage tank** provides a reserve, the pressure in an indirect-feed system can be increased. Bacterial control is easier in a direct-feed system.
- A tight-fitting lid and conical or bowl-shaped base; drain from its lowest point; venting through a 0.45 micrometer hydrophobic air filter.



2025 UpToDate, Inc. Topic 1896 Version 31.0

7- Distribution system components..


- The joints not to create any **rough areas** that can act as a focus for **bacterial proliferation**.
- The **number of fittings** in the circuit and abrupt changes in the direction of flow should be **minimized** as much as possible to prevent pockets of **stagnant flow** in the circuit.
- All surfaces within the distribution loop should be as smooth as possible .



Butt Fusion Pipe



A



B

Water Treatment Systems

1. Booster pump
2. Ultrafilters, particle filters, multimedia filters, and depth filters
3. Softeners
4. Activated carbon beds or filters
5. Reverse osmosis
6. Deionization
7. Distribution system components
- 8. Disinfection**

8- Disinfection

The risk of microbial recontamination in the distribution system is high due to **no chlorine and monochloramine** in the treatment water.

Frequency of disinfection: at a minimum should be performed **quarterly** even the **action level have not been breached**.

Disinfection:

1. **Chemicals:** sodium hypochlorite (bleach) or mixture of peracetic acid and hydrogen peroxide (compatible with the PVC).
2. Disinfection is more effective if pipes are first treated with a **descaling agent**.
3. **Ozonated water** or **hot water** is preferred over chemical disinfectants: lengthy period of rinsing required to remove residual disinfectant, no chemical residues (hot water) or residues with a short half-life (ozone), can be daily or weekly.

Water Monitoring systems



Water system monitoring 1

Water treatment system component	Sampling site	Test	Acceptable level	Recommended frequency of measurement
Feed water	Sampling port in the feed line to the water treatment system	Chemical analysis	Variable depending	<u>Every 3 months</u> or when the water provider changes the composition of the feed water or when the source of that water changes
Multimedia filters and/or beds	Pressure gauge before and after the filters and/or beds	Inflow versus outflow pressure difference and comparison with baseline pressure difference	>10 mmHg above baseline	<u>At the end of every day</u>
Softener	at the outflow of the softener	Hardness	<1 GPG (17 PPM)	<u>At the end of every day</u>

Water system monitoring 2

Water treatment system component	Sampling site	Test	Acceptable level	Recommended frequency of measurement
Carbon beds and filters	between the first and second carbon bed	Total chlorine Monochloramine	≤ 0.1 PPM	<u>At the start of the day and every 4 hours</u>
Reverse osmosis	N/A	Percent rejection	Based on water analysis, usually $>90\%$	<u>Continuously</u> (measured by reverse osmosis system).
Deionizer	N/A	DI resistivity	>1 megaohm•cm	<u>Continuously</u> (measured by deionizer)
Ultraviolet irradiator	N/A	Ultraviolet energy output	>30 mW•s/cm (>16 mW•s/cm if UV calibrated meter)	<u>Daily</u> , or in accordance with the manufacturer's recommendations if otherwise indicated

Water system monitoring 3

Water treatment system component	Sampling site	Test	Acceptable level	Recommended frequency of measurement
Dialysis water	Sampling port prior to the distribution loop	Chemical analysis	Variable, depending	<u>Every 3 months</u> or when: The quality of the feed water demonstrates changes in the measured chemical contaminants. The reverse osmosis membranes are changed. The reverse osmosis percent rejection drops below 90%.
	in the return line of the distribution loop in the distribution loop outlet to each treatment station	Colony count or total viable microbial count	<100 CFU/mL with an action level of 50 CFU/mL	From the return line of the distribution loop: <u>monthly</u> From the outlet of each treatment station: <u>annually</u>
	in the return line of the distribution loop , in the distribution loop outlet to each treatment station	Endotoxin	<0.25 EU/mL with an action level of 0.125 EU/mL	From the return line of the distribution loop: <u>monthly</u> From the outlet of each treatment station: <u>annually</u>



Review

Towards zero liquid discharge in hemodialysis. Possible issues

Faissal Tarrass*, Omar Benjelloun, Meryem Benjelloun

Center of Hemodialysis 2 Mars, Casablanca, Morocco

ARTICLE INFO

Article history:

Received 28 October 2020

Accepted 20 December 2020

ABSTRACT

Scarcity of water and energy, and legal requirements for discharge of waste and wastewater are forcing hemodialysis facilities to change their approach to a more integrated concept of connecting the residual output (in terms of waste, wastewater and energy loss) to the input (in terms of water and energy). Zero liquid discharge is an expanding water treat-

- **Every drop of water conserved reduces energy consumption and associated carbon emissions.**
- **Water use means energy use**, therefore conserving water means conserving energy.
- The carbon footprint of kidney care is high and estimated to be **27 million tonnes of carbon dioxide** per year worldwide.

NEFROLOGIA 2021 ;4 1(6):620–624



Advanced hemodialysis equipment for more eco-friendly dialysis

Adelheid Gaulty¹  · Nicole Fleck² · Fatih Kircelli¹

Received: 18 February 2021 / Accepted: 21 August 2021 / Published online: 4 September 2021

© The Author(s) 2021

Abstract

Healthcare in general and dialysis care in particular are contributing to resource consumption and, thus, have a notable environmental footprint. Dialysis is a life-saving therapy but it entails the use of a broad range of consumables generating waste, and consumption of water and energy for the dialysis process. Various stakeholders in the healthcare sector are called upon to develop and to take measures to save resources and to make healthcare and dialysis more sustainable. Among these stakeholders are manufacturers of dialysis equipment and water purification systems. Dialysis equipment and consumables, together with care processes need to be advanced to reduce waste generation, enhance recyclability, optimize water purification efficiency and water use. Joint efforts should thus pave the way to enable delivering green dialysis and to contribute to environmentally sustainable health care.

- It can be anticipated that with a dialysate flow of **500 ml/min** for 4 h not only the **120 L** water for **dialysis fluid** are required,
- But **up to 480 L** of water, taking also reject water into account.
- Additional water is used during the **preparation** and **disinfection** cycle adding to the water need by the actual treatment.

International Urology and Nephrology (2022) 54:1059–1065

Water implications in dialysis therapy, threats and opportunities to reduce water consumption: a call for the planet



OPEN

Mohamed Ben Hmida^{1,2}, Tahar Mechichi³, Giorgina Barbara Piccoli⁴ and Mohamed Ksibi⁵

¹Department of Nephrology, Hédi Chaker University Hospital, Sfax, Tunisia; ²Research Laboratory of Renal Pathology LR19ES11, Faculty of Medicine, University of Sfax, Sfax, Tunisia; ³Laboratory of Biochemistry and Enzymatic Engineering of Lipases, National Engineering School of Sfax, University of Sfax, Sfax, Tunisia; ⁴Department of Nephrology, Centre Hospitalier Le Mans, Le Mans, France; and ⁵Laboratory of Environmental Engineering and Eco-technology, National Engineering School of Sfax (LGEET-ENIS), University of Sfax, Sfax, Tunisia

Water is a dwindling natural resource, and potable water is wrongly considered an unlimited resource. Dialysis, particularly hemodialysis, is a water-hungry treatment that

Because of climate change, access to potable water is, and will be, a big challenge,¹ posing unprecedented threats to human health.² Traditionally, humans use water

Kidney International (2023) 104, 46–52; <https://doi.org/10.1016/j.kint.2023.04.008>

Hemodialysis: a water-hungry treatment

- The global annual **water use of hemodialysis** is approximately **265 million m³/yr.**
- “**3R**” concept:
 - **Reduce** (dialysis need, dialysate flow, optimize RO)
 - **Reuse** (reuse wastewater as potable water),
 - **Recycle** (dialysis effluents for agriculture and aquaponics use).
- Water consumption in hemodialysis is approximately **0.5 m³/session** based on the that two-thirds of this water is RO reject water discharged into the drain.
- Delaying dialysis initiation: for every patient-month of dialysis delay, the amount of water spared is approximately 6000 L (12 sessions x 500 L).

Kidney International (2023) 104, 46–52; <https://doi.org/10.1016/j.kint.2023.04.008>



Hemodialysis water reuse within a circular economy approach. What can we add to current knowledge? A point of view

Faissal Tarrass¹ · Meryem Benjelloun¹ · Giorgina Barbara Piccoli²

Received: 13 March 2024 / Accepted: 14 May 2024 / Published online: 3 June 2024
© The Author(s) 2024

Abstract

The ongoing climate change and the ecological challenges call for sustainable medicine and, in our field, sustainable kidney care. Dialysis is life-saving and resource-consuming, and high water consumption is one of the main concerns. Circular water economy, meaning reuse and recycling of water, and recovering resources can help reducing emissions and enhancing resilience to climate change. Several actions are possible including reusing reverse osmosis reject water, employable for gardening, aquaponics or even simply for toilet flushing, or in sterilization settings, reusing spent dialysate, at least for toilet flushing, but with wider use if microbiologically purified, recovering thermal energy from spent dialysate, that can probably be done with simple devices, or using phosphate-rich spent dialysate for producing fertilizers, namely struvite. All these options may be economically sound, and all help reducing the final dialysis carbon footprint. There is room for open-minded innovative approaches to improve water-related sustainability in hemodialysis, ultimately reducing ecological impact and increasing availability.

- Thermal energy is lost in dialysis units worldwide, could heat approximately 140,000 homes in a european country for one year.
- Spent dialysis effluent, which has a high saline content and is a potential source of infectious contamination, is systematically discharged.
- Since an excessive salt level can have a negative effect on plant growth, mixing reverse osmosis reject water with water from other sources (rain, wells) allows us to dilute salts to a safe level for use in farming.
- Reverse osmosis and nanofiltration, are effective in reducing not only salt content but also an array of micropollutants and pathogens, including viruses, drugs and drug metabolites.

Journal of Nephrology (2024) 37:1801–1805

Case cost studies of reject reverse osmosis water reuse in hemodialysis

	Canterbury Hospital NHS—UK [9]	Countess of Chester Hospital NHS—UK [10]	Lister Hospital NHS—UK [11]	Maidland Hospital Ireland [12]	Sultan Abdul Halim Hospital Malaysia [13]
Reuse option	Toilet flushing	Toilet flushing, Laundry and other uses	Hot water	Toilet flushing	Acquaponics and horti- culture
Volume of saved water	800 Liters / hour	1460 M ³ /year	3140 m ³ / year	5240 m ³ / year	12,000 liters/day
Implementation costs	19,000 \$ (15,000 £)	14,000 \$ (11,030 £)	7600 \$ (6000 £)	13,100 \$ (12,000 €)	\$1047 (MYR 5000)
Financial savings (USD/year)	9530 \$ (7500 £)	3990 \$ (3143 £)	8000 \$ (6300 £)	13,600 \$ (12,500 €)	\$524 (MYR 2500)
Payback period (months)	24	42	12	12	24

Take home messages

1. The major components of water treatment system include: multimedia or sediment filter, water softener, carbon filter and RO.
2. **Multimedia filter** removes particles > 20 micron
3. **Water softener** removes calcium and magnesium.
4. **Carbon filter** removes chlorine and monochloramine.
5. **Reverse osmosis** removes metal ions, salts, inorganic, organic chemicals and microbial elements.
6. Materials of distribution system should be from **PEX, PVC or stainless steel**.

Take home messages

7. Bacterial analysis from outlet of each dialysis machine annually.
8. Feed and distribution water chemical analysis every 3 months.
9. Bacterial analysis of dialysis water every months.
10. Sediment filter, softener monitoring at the end of every day.
11. Carbon filter at the beginning and every 4 hours.
12. Water consumption in hemodialysis is approximately 0.5 m³/session.
13. The carbon footprint of kidney care is high and estimated to be 27 million tonnes of carbon dioxide per year worldwide.